



**FRESENIUS
KABI**

***** **CERTIFICATE OF ANALYSIS** *****

Article: CompoStop Flex2F PLT processing and leuko red.syst
 Article number: FD51600
 Batch number: 41WI30GA00 Quantity: 15 320
 Manufacturing date: 29 September 2024 Expiry date: 2026-08

Test	Limits	Results
Configuration	Conform product specification	Conform
Physical control	In accordance with specifications	Conform
Bacterial Endotoxines	< 20 EU/device	Conform
Sterility	Sterile (Ph.Eur.)	Conform
Used Filterbatch	24PA015	

<u>Batch released by Quality Assurance.</u>	<u>Manufactured by:</u>
Date: 16.10.2024	Fresenius Kabi AG
Name: Klímová Vladimíra	Plant:
Signature: <i>V. Klímová</i>	Fresenius Kabi Horatev CZ s.r.o.
	Hořátev 104, 289 13
	Czech Republic