



**FRESENIUS  
KABI**

\*\*\*\*\* **CERTIFICATE OF ANALYSIS** \*\*\*\*\*

Article: CompoStop Flex2F PLT processing and leuko red.sys  
Article number: FD51600  
Batch number: 41WL21GA00 Quantity: 15 480  
Manufacturing date: 18 November 2024 Expiry date: 2026-10

<b>Test</b>	<b>Limits</b>	<b>Results</b>
Configuration	Conform product specification	Conform
Physical control	In accordance with specifications	Conform
Bacterial Endotoxines	< 20 EU/device	Conform
Sterility	Sterile (Ph.Eur.)	Conform
Used Filterbatch	24PA021,24PA022	

<u>Batch released by Quality Assurance.</u>	<u>Manufactured by:</u>
Date: 06.12.2024	Fresenius Kabi AG
Name: Tykva František	Plant: Fresenius Kabi Horatev CZ s.r.o.
Signature: 	Hořátev 104, 289 13 Czech Republic